

## SmartVision2

### Release note – software version n° SV2-7602

#### Software Update Process

An Internet connection is required to download and install SmartVision2 updates. To update your device, go to “Settings”, select "About Phone" and then "System Updates". Select "Download and install" at the bottom of the list and press the "OK" button to start the download.

Depending on your connection, downloading the update may take some time. You can view its progress via the update notification in the notification shade. When the download is complete, select "Install" from the system update screen. It may take a few minutes to complete the process before you could use the product (vocalized feedback).

#### System improvements:

- **Overview:**
  - Bugs fix and performance improvement
  - Translation and vocalization improvement
  - New languages: Croatian and Latvian
- **System:**
  - Copy / Cut / Paste will not store anymore the last vocal recognition result (Keypad shortcut \*+3 (press and hold \* key and then press 3), \*+6 and \*+9).
  - Voice dictation fixed bug: no more double text insertion within input boxes during voice dictation (due to a Google update)
- **Accessibility:**
  - Improvement of vocalization for both PIN code and wake-up speeches when the "Text-to-speech Multi-Language" mode is enabled.
  - Reduced verbosity when using the shortcut #+1 (press and hold # key and then press 3) to get phone status.
  - Improvement of the Magnitext display in landscape mode

#### Application Improvements:

- **Assistant**
  - Activating Google Assistant function with the voice recognition button is now faster and more precise.
- **Calendar**
  - New "Save" button at the end of the list to simplify the creation of an event.
  - New easier way to delete an event from the list: press on “Delete” key
  - The most recent event is now displayed first when the application is opened

- The current date and hour are now displayed first when you create a new event
- Fixed bug when deleting an event created from the Google Calendar application
- **Contacts**
  - The confirmation pop-up message for creating a new contact is now vocalized.
- **Banknote recognizer**
  - The "Help" item from Menu (press Menu key) is now available to allow the user access to the Banknote Recognizer chapter of the User Manual
  - Flash LED is automatically activated to improve the detection accuracy of banknotes.
- **Sound Recorder**
  - Voice memo are now sorted from newest to oldest.
- **Email**
  - New option to download automatically the full content of an Email (previous discussions and attachments).
  - When you delete an Email with the "Delete" key, the focus is now placed correctly on the next Email.
- **Kapten**
  - Titles of each screen are now announced systematically.
  - The warning message about KAPTEN use is now announced only one time at the first activation
- **Magnifier**
  - Magnifier application is redesigned using a new image processing method resulting significant quality improvements as follow:
    - Fluid image display
    - better stabilization
    - enhanced contours
    - increased zoom level
  - The favorite settings of the user are now saved for the next use. No more need to set them again.
  - New feature: choose to start yourself the camera focus whenever you are ready by pressing the OK button.
- **Audio Book**
  - This reader application is redesigned using a new engine resulting multiple improvements and bug fixes.
  - New easier way to continue the reading of a current book: All books are now resumed at their last position when reopened.

- New feature to remove easily a book:  
Select a book from the bookshelf, make a long press on “OK” key and select “Remove”. A confirmation pop-up will confirm the suppression.
- Improved reading for different formats: Daisy 2.02, Daisy 3.0, Epub, Epub3, Doc, Docx PDF, Txt et html.
- Quick opening of large books or files
- Issue correction on “Resume” feature for some specific books or files
  
- **Message**
  - New feature for MMS:  
Possibility to send several pictures in the same MMS message
  - More than 50 new emojis available
  - Bug fix related to "Search" function when user have lots of thread conversations
  
- **OCR**
  - Improved online detection mode
  - Improved detection of accented characters and symbols
  
- **Phone**
  - New: possibility to call directly the last phone number with a long press on “Pick-up” key.
  - New: possibility to remove an entry from the call history list with “Delete” key.
  - The “Phone application has stopped” issue has been fixed.
  
- **QuickSupport**
  - The application is updated to assure the compatibility with the last version of TeamViewer